
Oncomine Precision Assay on Genexus NGS System

A next-generation solution for genomic profiling

Empower **your lab** to deliver NGS genomic profiles with the **speed** and **simplicity** of immunohistochemistry. Maximize your ability to deliver results even from small samples and low-level variants.

Minimum sample requirement



Maximum results obtained



Single-day sample-to-report workflow* 2 touch points and only 10 min of hands-on time






Curated pan-cancer content

- Mutations, CNVs, and fusion variant types across 50 key genes
- Tumor suppressors, drivers, and resistance variants



Tissue and plasma samples

- One test, one workflow, multiple sample types
- Maximizes the number of tumors that can be profiled



Molecular tagging

- Enhanced low-level variant detection
- Key for liquid biopsy testing



FusionSync™ Detection Technology

- Sensitive and specific—targeted isoform designs
- Novel fusion detection

Find out more at www.oncomine.com/genexus-oncology



^{*} Specimen-to-report workflow will be available after the Genexus Purification System and integrated reporting capabilities are added in 2020.